

Product Datasheet

AMPK alpha 1 Antibody orb1517493

Description	AMPK alpha 1 Antibody
Species/Host	Rabbit
Reactivity	Bovine, Human, Rat
Conjugation	Agarose
Tested Applications	IP
Immunogen	Between 350 and 400
Target	AMPK alpha 1
Form/Appearance	Whole IgG
Preservatives	Phosphate Buffered Saline (PBS) containing 0.09% Sodium Azide
Storage	2 - 8°C
Note	For research use only.
Clonality	Polyclonal
Uniprot ID	Q13131
NCBI	BAA36547.1
Dilution Range	IP - 15 - 25 µl of gel slurry per 0.1 to 1 mg of protein lysate or extract
Expiration Date	12 months from date of receipt.