

Product Datasheet

MAPT Antibody orb1517462

Description	MAPT Monoclonal Antibody
Reactivity	Monkey, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Human MAPT protein
Target	MAPT
Form/Appearance	Liquid
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only.
Isotype	hIgG4(S228P)
Clonality	Monoclonal
Clone ID	8C2
Uniprot ID	P10636
Expiration Date	12 months from date of receipt.