

Product Datasheet

Envelope protein p54 Antibody orb1517349

Description	Ba71V-126 Monoclonal Antibody
Species/Host	Mouse
Reactivity	Virus
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	Recombinant Envelope protein p54 protein (54-183AA)
Target	Envelope protein p54
Form/Appearance	Liquid
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only.
Clonality	Monoclonal
Clone ID	2F7H5
Uniprot ID	Q65194
Expiration Date	12 months from date of receipt.