

## Product Datasheet

### African swine fever virus Phosphoprotein p30 Antibody orb1517339

|                            |                                                                            |
|----------------------------|----------------------------------------------------------------------------|
| <b>Description</b>         | African swine fever virus Phosphoprotein p30(Ba71V-93) Monoclonal Antibody |
| <b>Species/Host</b>        | Mouse                                                                      |
| <b>Reactivity</b>          | Virus                                                                      |
| <b>Conjugation</b>         | Unconjugated                                                               |
| <b>Tested Applications</b> | ELISA                                                                      |
| <b>Immunogen</b>           | Recombinant Phosphoprotein p30 protein (1-204AA)                           |
| <b>Target</b>              | African swine fever virus Phosphoprotein p30                               |
| <b>Form/Appearance</b>     | Liquid                                                                     |
| <b>Preservatives</b>       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol                |
| <b>Storage</b>             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.              |
| <b>Note</b>                | For research use only.                                                     |
| <b>Isotype</b>             | IgG1                                                                       |
| <b>Clonality</b>           | Monoclonal                                                                 |
| <b>Clone ID</b>            | 5F2F3                                                                      |
| <b>Uniprot ID</b>          | <b>P34204</b>                                                              |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                            |