

## Product Datasheet

### African swine fever virus Phosphoprotein p30 Antibody orb1517338

|                     |                                                                            |
|---------------------|----------------------------------------------------------------------------|
| Description         | African swine fever virus Phosphoprotein p30(Ba71V-93) Monoclonal Antibody |
| Species/Host        | Mouse                                                                      |
| Reactivity          | Virus                                                                      |
| Conjugation         | Unconjugated                                                               |
| Tested Applications | ELISA                                                                      |
| Immunogen           | Recombinant Phosphoprotein p30 protein (1-204AA)                           |
| Target              | African swine fever virus Phosphoprotein p30                               |
| Form/Appearance     | Liquid                                                                     |
| Preservatives       | Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol                |
| Storage             | Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.              |
| Note                | For research use only.                                                     |
| Isotype             | IgG1                                                                       |
| Clonality           | Monoclonal                                                                 |
| Clone ID            | 2E1F3                                                                      |
| Uniprot ID          | <b>P34204</b>                                                              |
| Expiration Date     | 12 months from date of receipt.                                            |