

Product Datasheet

M Monoclonal Antibody orb1517337

Description	M Monoclonal Antibody
Species/Host	Mouse
Reactivity	Rabies virus
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Rabies virus Matrix protein (1-202AA)
Target	M
Form/Appearance	Liquid
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only.
Isotype	IgG2b
Clonality	Monoclonal
Clone ID	1F9H6
Uniprot ID	P15200
Expiration Date	12 months from date of receipt.