

Product Datasheet

III Antibody orb1517336

Description	Enterobacteria phage M13 III Monoclonal Antibody
Species/Host	Mouse
Reactivity	Bacteria
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	Recombinant Enterobacteria phage M13 Attachment protein (19-424AA)
Target	III
Form/Appearance	Liquid
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only.
Isotype	IgG1
Clonality	Monoclonal
Clone ID	4B1D7
Uniprot ID	P69168
Expiration Date	12 months from date of receipt.