

Product Datasheet

CD84 polyclonal antibody?KU? orb1496109

Description	CD84 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Tested Applications	WB
Immunogen	Recombinant protein of human CD84.
Concentration	1mg/ml
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Polyclonal
MW	~ 60 kDa
Uniprot ID	Q9UIB8
Dilution Range	WB: 1:1000~1:2000
Expiration Date	12 months from date of receipt.