

Product Datasheet

FAM134B polyclonal antibody orb1496092

Description	FAM134B polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Tested Applications	ELISA, ICC, IF, IHC
Immunogen	KLH conjugated synthetic peptide derived from human FAM134B: 210-310/497
Concentration	1mg/ml
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Polyclonal
MW	~ 55 kDa
Uniprot ID	Q9H6L5
Dilution Range	ELISA: 1:5000-1:10000 IHC: 1:100-1:500 ICC: 1:100-1:500 IF: 1:100-1:500
Expiration Date	12 months from date of receipt.