

Product Datasheet

GSK3 β (Ab-9) polyclonal antibody (KU) orb1496046

Description	GSK3 β (Ab-9) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Tested Applications	IF, WB
Immunogen	Synthetic peptide, corresponding to Human GSK3 β .
Concentration	1mg/ml
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Polyclonal
MW	~ 46 kDa
Uniprot ID	P49841
Dilution Range	WB: 1:1000~1:2000
Expiration Date	12 months from date of receipt.