

Product Datasheet

ATG4D (G30) polyclonal antibody orb1496037

Description	ATG4D (G30) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Tested Applications	WB
Immunogen	Synthetic peptide, corresponding to Human ATG4D.
Concentration	1mg/ml
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Polyclonal
MW	~ 52 kDa
Uniprot ID	Q86TL0
Dilution Range	WB: 1:1000~1:2000
Expiration Date	12 months from date of receipt.