

Product Datasheet

LKB1(Phospho-Ser428) polyclonal antibody?KU? orb1496034

Description	LKB1(Phospho-Ser428) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Tested Applications	WB
Immunogen	Synthetic phosphopeptide derived from human LKB1 around the phosphorylation site of Serine 428.
Concentration	1mg/ml
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Polyclonal
MW	~ 55 kDa
Uniprot ID	Q15831
Dilution Range	WB: 1:1000~1:2000
Expiration Date	12 months from date of receipt.