

Product Datasheet

ATG9B (G883) polyclonal antibody?KU? orb1496024

Description	ATG9B (G883) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Tested Applications	IP, WB
Immunogen	Synthetic peptide, corresponding to Human ATG9B.
Concentration	1mg/ml
Preservatives	Rabbit IgG, 1mg/ml in PBS with 0.02% sodium azide, 50% glycerol, pH 7.2.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Polyclonal
MW	~ 95 kDa
Uniprot ID	Q674R7
Dilution Range	WB: 1:5000~1:10000 IP: 1:50~1:200
Expiration Date	12 months from date of receipt.