

Product Datasheet

Histone Deacetylase 6 (phospho-S22) polyclonal antibody orb1495958

Description	Histone Deacetylase 6 (phospho-S22) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human
Tested Applications	IF, IHC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human Histone Deacetylase 6 (pS22). The exact sequence is proprietary.
Concentration	1mg/ml
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Polyclonal
MW	~ 160 kDa
Uniprot ID	Q9UBN7
Dilution Range	WB (1/500 - 1/1000), IH (1/100 - 1/200), IF/IC (1/100 - 1/500)
Expiration Date	12 months from date of receipt.