

## Product Datasheet

### USP51 polyclonal antibody orb1495883

<b>Description</b>	USP51 polyclonal antibody
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human
<b>Tested Applications</b>	IHC, WB
<b>Immunogen</b>	KLH-conjugated synthetic peptide encompassing a sequence within the center region of human USP51. The exact sequence is proprietary.
<b>Concentration</b>	1mg/ml
<b>Preservatives</b>	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
<b>Note</b>	For research use only.
<b>Clonality</b>	Polyclonal
<b>MW</b>	~ 80 kDa
<b>Uniprot ID</b>	<b>Q70EK9</b>
<b>Dilution Range</b>	WB (1/500 - 1/1000), IHC (1/50 - 1/200)
<b>Expiration Date</b>	12 months from date of receipt.