

Product Datasheet

WIPI2 polyclonal antibody orb1495831

Description	WIPI2 polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the C-term region of human WIPI2. The exact sequence is proprietary.
Concentration	1mg/ml
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Polyclonal
MW	~ 50 kDa
Uniprot ID	Q9Y4P8
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.