

Product Datasheet

HMGB1 (AcK12) polyclonal antibody orb1495767

Description	HMGB1 (AcK12) polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Tested Applications	WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within the N-term region of human HMGB1 (AcK12). The exact sequence is proprietary.
Concentration	1mg/ml
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Polyclonal
MW	~ 90 kDa
Uniprot ID	P09429
Dilution Range	WB (1/500 - 1/1000)
Expiration Date	12 months from date of receipt.