

Product Datasheet

CD183 Armenian Hamster monoclonal antibody orb1495658

Description	CD183 Armenian Hamster monoclonal antibody
Species/Host	Hamster
Reactivity	Mouse
Tested Applications	FC, IF
Immunogen	Mouse N-terminus of CXCR3
Concentration	1mg/ml
Preservatives	Armenian Hamster IgG. Liquid in PBS, pH 7.3, and 0.02% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Monoclonal
MW	~ 41 kDa
Uniprot ID	O88410
Dilution Range	IF (1/50 - 1/200)
Expiration Date	12 months from date of receipt.