

Product Datasheet

Cleaved PARP1 monoclonal antibody orb1495592

Description	Cleaved PARP1 monoclonal antibody
Species/Host	Mouse
Reactivity	Human
Tested Applications	IHC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of human Cleaved PARP1. The exact sequence is proprietary.
Concentration	1mg/ml
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Monoclonal
MW	~ 89 kDa
Uniprot ID	P09874
Dilution Range	WB (1/2000 - 1/5000), IHC (1/200 - 1/500)
Expiration Date	12 months from date of receipt.