

Product Datasheet

NF-kappaB p65 monoclonal antibody orb1495533

Description	NF-kappaB p65 monoclonal antibody
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Tested Applications	IF, IHC, IP, WB
Immunogen	Recombinant protein corresponding to human NF-kappaB p65.
Concentration	1mg/ml
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Monoclonal
MW	~ 65 kDa
Uniprot ID	Q04206
Dilution Range	WB (1/500 - 1/1000), IHC (1/200 - 1/500), IF/ICC (1/100 - 1/200), IP (1/100 - 1/200)
Expiration Date	12 months from date of receipt.