

Product Datasheet

Ghrelin Receptor polyclonal antibody orb1495510

Description	Ghrelin Receptor polyclonal antibody
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Tested Applications	IHC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence of human Ghrelin Receptor. The exact sequence is proprietary.
Concentration	1mg/ml
Preservatives	Liquid in 0.42% Potassium phosphate, 0.87% Sodium chloride, pH 7.3, 30% glycerol, and 0.01% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Polyclonal
MW	~ 41 kDa
Uniprot ID	Q92847
Dilution Range	WB (1/1000 - 1/2000), IHC (1/100 - 1/200)
Expiration Date	12 months from date of receipt.