

Product Datasheet

Histone Deacetylase 6 monoclonal antibody orb1495468

Description	Histone Deacetylase 6 monoclonal antibody
Species/Host	Mouse
Reactivity	Human
Tested Applications	IHC
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human Histone Deacetylase 6. The exact sequence is proprietary.
Concentration	1mg/ml
Preservatives	Mouse IgG1. Liquid in PBS containing 50% glycerol, 0.2% rAlbumin and 0.01% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Monoclonal
MW	~ 131 kDa
Uniprot ID	Q9UBN7
Dilution Range	IHC (1/100 - 1/300)
Expiration Date	12 months from date of receipt.