

Product Datasheet

Actin Muscle Specific monoclonal antibody orb1495458

Description	Actin Muscle Specific monoclonal antibody
Species/Host	Mouse
Reactivity	Human
Tested Applications	IHC
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human Actin Muscle Specific. The exact sequence is proprietary.
Concentration	1mg/ml
Preservatives	Mouse IgG. Liquid in PBS containing 50% glycerol, 0.2% rAlbumin and 0.01% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Monoclonal
MW	~ 42 kDa
Uniprot ID	P68032 , P63267 , P68133
Dilution Range	IHC (1/100 - 1/300)
Expiration Date	12 months from date of receipt.