

Product Datasheet

Desmin monoclonal antibody orb1495428

Description	Desmin monoclonal antibody
Species/Host	Mouse
Reactivity	Human
Tested Applications	IHC
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human Desmin. The exact sequence is proprietary.
Concentration	1mg/ml
Preservatives	Mouse IgG2b. Liquid in PBS containing 50% glycerol, 0.2% rAlbumin and 0.01% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Monoclonal
MW	~ 53 kDa
Uniprot ID	P17661
Dilution Range	IHC (1/100 - 1/300)
Expiration Date	12 months from date of receipt.