

Product Datasheet

Claudin 4 monoclonal antibody orb1495337

Description	Claudin 4 monoclonal antibody
Species/Host	Mouse
Reactivity	Human
Tested Applications	IHC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human Claudin 4. The exact sequence is proprietary.
Concentration	1mg/ml
Preservatives	Mouse IgG2b. Liquid in PBS containing 50% glycerol, 0.2% rAlbumin and 0.01% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Monoclonal
MW	~ 23 kDa
Uniprot ID	O14493
Dilution Range	WB (1/500 - 1/1000), IHC (1/100 - 1/300)
Expiration Date	12 months from date of receipt.