

Product Datasheet

Actin Sarcomeric Muscle monoclonal antibody orb1495274

Description	Actin Sarcomeric Muscle monoclonal antibody
Species/Host	Mouse
Reactivity	Human
Tested Applications	IHC, WB
Immunogen	KLH-conjugated synthetic peptide encompassing a sequence within human Actin Sarcomeric Muscle. The exact sequence is proprietary.
Concentration	1mg/ml
Preservatives	Mouse IgG1. Liquid in PBS containing 50% glycerol, 0.2% rAlbumin and 0.01% sodium azide.
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Note	For research use only.
Clonality	Monoclonal
MW	~ 42 kDa
Uniprot ID	P68032 , P68133
Dilution Range	WB (1/500 - 1/1000), IHC (1/100 - 1/300)
Expiration Date	12 months from date of receipt.