

Product Datasheet

CD312 Recombinant Protein orb1495113

Description	CD312 Recombinant Protein
Species/Host	E. coli
Preservatives	PBS, 4M Urea, pH 7.4
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze-thaw cycles.
Tag	His-tag
Note	For research use only.
Protein Sequence	QDSRGCARWCPQDSSCVNATACRCNPGFSSFSEIITTPMETCDDINECAT LSKVSCGKFSDCWNTGSGYDCVCSPGYEPVSGAKTFKNESENTCQDVDEC QQNPRLCKSYGTCVNTLGSYTCQCLPGFKLPEDPKLCTDVNECTSGQNP CHSSTHCLNNVGSYQCRCPGWQPIPGSPNGPNNTVCEDVDECSSGQHQC DSSTVCFNTVGSYSCRCRPGWKPRHGIPNNQKDTVCEDMTFSTWTPPPGV HSQTLRFFDKVQDLGRDYKPLANNTIQSILQALDELLEAPGDLETLP LQQHCVASHLLDGLDVLRLSKNLSNGLLNFSYPAGTELSLEVQKQVDR SVTLRQNQAVMQLDWNQAQKSGDPGPSVVGLVSIPGMGKLLAEAPLVLEP EKQMLLHETHQGLLQDGSPILLSDVISAFLSNNDTQNLSSPVTFTFSHRS VIPRQKVLGVFWEHGQNGCGHWATTGCSTIGTRDTSTICRCTHLSSFAVL MAHYDVQEEDPVLTVIT
MW	~55kDa
Uniprot ID	Q9UHX3
Solubility (25°C)	PBS, 4M Urea, PH7.4
Expiration Date	6 months from date of receipt.