

## Product Datasheet

### Recombinant Human CD120A/TNFRSF1A Protein, His Tag orb1553328

<b>Description</b>	Recombinant Human CD120A/TNFRSF1A Protein, His Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Human
<b>Endotoxins</b>	
<b>Target</b>	TNFRSF1A
<b>Note</b>	For research use only.
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human TNFRSF1A(Met1-Thr211) was fused with the C-terminal His Tag and expressed in Mammalian cells.
<b>Uniprot ID</b>	<b>P19438</b>
<b>Expiration Date</b>	6 months from date of receipt.