

## Product Datasheet

### Recombinant Human CD121A/IL1R1 Protein, Fc Tag orb1553326

<b>Description</b>	Recombinant Human CD121A/IL1R1 Protein, Fc Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Human
<b>Endotoxins</b>	
<b>Target</b>	IL1R1
<b>Note</b>	For research use only.
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human IL1R1(Met1-Thr332) was fused with the C-terminal Fc Tag and expressed in Mammalian cells.
<b>Uniprot ID</b>	<a href="#">P14778</a>
<b>Expiration Date</b>	6 months from date of receipt.