

Product Datasheet

Recombinant Human TNF Protein, His Tag orb1553128

Description	Recombinant Human TNF Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	
Target	TNF
Note	For research use only.
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human TNF(Val77-Leu233) was fused with the N-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	P01375
Expiration Date	6 months from date of receipt.