

## Product Datasheet

### Recombinant Staphylococcus aureus Protein A/spa Protein, His Tag orb1553022

<b>Description</b>	Recombinant Staphylococcus aureus Protein A/spa Protein, His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Bacteria
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	spa
<b>Note</b>	For research use only.
<b>Purity</b>	?90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the Staphylococcus aureus spa(Met1-Lys327) was fused with the C-terminal His Tag and expressed in E. coli.
<b>Uniprot ID</b>	<b>P38507</b>
<b>Expiration Date</b>	6 months from date of receipt.