

## Product Datasheet

### Recombinant *Streptococcus pyogenes* serotype M1 cas9 Protein, His tag orb1553001

<b>Description</b>	Recombinant <i>Streptococcus pyogenes</i> serotype M1 cas9 Protein, His tag
<b>Species/Host</b>	<i>E. coli</i>
<b>Reactivity</b>	Bacteria
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	cas9
<b>Note</b>	For research use only.
<b>Purity</b>	≥90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the <i>Streptococcus pyogenes</i> serotype M1 cas9(Met1-Asp1368) was fused with the N-terminal His Tag . and expressed in <i>E. coli</i> .
<b>Uniprot ID</b>	<a href="#">Q99ZW2</a>
<b>Expiration Date</b>	6 months from date of receipt.