

Product Datasheet

Recombinant Human ATG4C Protein, His Tag orb1552889

Description	Recombinant Human ATG4C Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	ATG4C
Note	For research use only.
Purity	?90% as determined by SDS-PAGE
Source	This protein contains the human ATG4C(Met1-Leu458) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	Q96DT6
Expiration Date	6 months from date of receipt.