

## Product Datasheet

### Recombinant Human BAI3/ADGRB3 Protein, His Tag orb1552863

Description	Recombinant Human BAI3/ADGRB3 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	
Target	ADGRB3
Note	For research use only.
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human ADGRB3(Ala25-Thr880) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	<a href="#">O60242</a>
Expiration Date	6 months from date of receipt.