

Product Datasheet

Recombinant Human Cadherin-3/CALR3 Protein, No Tag orb1552786

Description	Recombinant Human Cadherin-3/CALR3 Protein, No Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	CALR3
Note	For research use only.
Purity	?90% as determined by SDS-PAGE
Source	This protein contains the human CALR3(Thr20-Leu384) was fused without Tag and expressed in E. coli.
Uniprot ID	Q96L12
Expiration Date	6 months from date of receipt.