

Product Datasheet

Recombinant Human CDK7 Protein, His Tag orb1552711

Description	Recombinant Human CDK7 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	CDK7
Note	For research use only.
Purity	?90% as determined by SDS-PAGE
Source	This protein contains the human CDK7(Met1-Phe346) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	P50613
Expiration Date	6 months from date of receipt.