

Product Datasheet

Recombinant Human DNA Polymerase/POLD2 Protein, His Tag orb1552535

| | |
|------------------------|--|
| Description | Recombinant Human DNA Polymerase/POLD2 Protein, His Tag |
| Species/Host | E. coli |
| Reactivity | Human |
| Endotoxins | Please contact our customer service for more information. |
| Target | POLD2 |
| Note | For research use only. |
| Purity | ?90% as determined by SDS-PAGE |
| Source | This protein contains the human POLD2(Met1-Pro469) was fused with the C-terminal His Tag and expressed in E. coli. |
| Uniprot ID | P49005 |
| Expiration Date | 6 months from date of receipt. |