

Product Datasheet

Recombinant Human FE1/FEN1 Protein, His Tag orb1552422

Description	Recombinant Human FE1/FEN1 Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	FEN1
Note	For research use only.
Purity	?90% as determined by SDS-PAGE
Source	This protein contains the human FEN1(Met1-Lys380) was fused with the C-terminal His Tag and expressed in E. coli.
Uniprot ID	P39748
Expiration Date	6 months from date of receipt.