

## Product Datasheet

### Recombinant Human GABARAP Protein, His Tag & Fc Tag orb1552360

<b>Description</b>	Recombinant Human GABARAP Protein, His Tag & Fc Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	GABARAP
<b>Note</b>	For research use only.
<b>Purity</b>	?90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human GABARAP(Met1-Gln116) was fused with the N-terminal His Tag,C-terminal Fc Tag and expressed in E. coli.
<b>Uniprot ID</b>	<b>O95166</b>
<b>Expiration Date</b>	6 months from date of receipt.