

## Product Datasheet

### Recombinant Human GADD45A Protein, His Tag orb1552357

Description	Recombinant Human GADD45A Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	GADD45A
Note	For research use only.
Purity	?90% as determined by SDS-PAGE
Source	This protein contains the human GADD45A(Met1-Arg165) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	<b>P24522</b>
Expiration Date	6 months from date of receipt.