

Product Datasheet

Recombinant Human HDAC10 Protein, His Tag orb1552273

Description	Recombinant Human HDAC10 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	
Target	HDAC10
Note	For research use only.
Purity	?90% as determined by SDS-PAGE
Source	This protein contains the human HDAC10(Met1-Ala669) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	Q969S8
Expiration Date	6 months from date of receipt.