

## Product Datasheet

### Recombinant Human HMGB1 Protein, MARI Tag orb1552260

|                 |  |
|-----------------|--|
| Description     | Recombinant Human HMGB1 Protein, MARI Tag  |
| Species/Host    | E. coli  |
| Reactivity      | Human  |
| Endotoxins      | Please contact our customer service for more information.  |
| Target          | HMGB1  |
| Note            | For research use only.   |
| Purity          | ?90% as determined by SDS-PAGE   |
| Source          | This protein contains the human HMGB1(Pro92-Val176) was fused with the N-terminal MARI Tag and expressed in E. coli. |
| Uniprot ID      | <b>P09429</b>  |
| Expiration Date | 6 months from date of receipt.   |