

Product Datasheet

Recombinant Mouse Il-1 RI/Il1r1 Protein, Fc Tag orb1552198

Description	Recombinant Mouse Il-1 RI/Il1r1 Protein, Fc Tag
Species/Host	Mammal
Reactivity	Mouse
Endotoxins	
Target	Il1r1
Note	For research use only.
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Il1r1(Leu20-Lys338) was fused with the C-terminal Fc Tag and expressed in Mammalian cells.
Uniprot ID	P13504
Expiration Date	6 months from date of receipt.