

Product Datasheet

Recombinant Mouse Il-4 Protein, His Tag orb1552174

Description	Recombinant Mouse Il-4 Protein, His Tag
Species/Host	Mammal
Reactivity	Mouse
Endotoxins	
Target	Il4
Note	For research use only.
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the mouse Il4(His21-Ser140) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	P07750
Expiration Date	6 months from date of receipt.