

## Product Datasheet

### Recombinant Mouse Mrp8/S100a8 Protein, His Tag orb1552000

<b>Description</b>	Recombinant Mouse Mrp8/S100a8 Protein, His Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	S100a8
<b>Note</b>	For research use only.
<b>Purity</b>	?90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the mouse S100a8(Met1-Glu93) was fused with the C-terminal His Tag and expressed in E. coli.
<b>Uniprot ID</b>	<b>P05109</b>
<b>Expiration Date</b>	6 months from date of receipt.