

Product Datasheet

Recombinant *Streptomyces hygroscopicus* Phosphinothricin N-acetyltransferase/bar Protein, No Tag orb1551843

Description	Recombinant <i>Streptomyces hygroscopicus</i> Phosphinothricin N-acetyltransferase/bar Protein, No Tag
Species/Host	<i>E. coli</i>
Reactivity	Bacteria
Endotoxins	Please contact our customer service for more information.
Target	PPT N-acetyltransferase
Note	For research use only.
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the <i>Streptomyces hygroscopicus</i> bar(Met1-Ile183) was fused without Tag and expressed in <i>E. coli</i> .
Uniprot ID	P16426
Expiration Date	6 months from date of receipt.