

Product Datasheet

Recombinant Human PTH1R/PTHR1 Protein, His Tag orb1551767

Description	Recombinant Human PTH1R/PTHR1 Protein, His Tag
Species/Host	Mammal
Reactivity	Human
Endotoxins	
Target	PTH1R
Note	For research use only.
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human PTH1R(Tyr23-Met189) was fused with the C-terminal His Tag and expressed in Mammalian cells.
Uniprot ID	Q03431
Expiration Date	6 months from date of receipt.