

Product Datasheet

Recombinant Human Synuclein α /SNCA N-his Protein, His Tag orb1551574

Description	Recombinant Human Synuclein α /SNCA N-his Protein, His Tag
Species/Host	E. coli
Reactivity	Human
Endotoxins	Please contact our customer service for more information.
Target	SNCA
Note	For research use only.
Purity	≥90% as determined by SDS-PAGE
Source	This protein contains the human SNCA(Met1-Ala140) was fused with the N-terminal His Tag and expressed in E. coli.
Uniprot ID	P37840
Expiration Date	6 months from date of receipt.