

## Product Datasheet

# Recombinant Human Synuclein $\alpha$ /SNCA Protein, No Tag orb1551573

<b>Description</b>	Recombinant Human Synuclein $\alpha$ /SNCA Protein, No Tag
<b>Species/Host</b>	E. coli
<b>Reactivity</b>	Human
<b>Endotoxins</b>	Please contact our customer service for more information.
<b>Target</b>	SNCA
<b>Note</b>	For research use only.
<b>Purity</b>	?90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the human SNCA(Met1-Ala140) was fused without Tag and expressed in E. coli.
<b>Uniprot ID</b>	<b>P37840</b>
<b>Expiration Date</b>	6 months from date of receipt.