

## Product Datasheet

### Recombinant Mouse TGF-beta RII/Tgfbr2 Protein, Fc Tag orb1551557

<b>Description</b>	Recombinant Mouse TGF-beta RII/Tgfbr2 Protein, Fc Tag
<b>Species/Host</b>	Mammal
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	
<b>Target</b>	Tgfbr2
<b>Note</b>	For research use only.
<b>Purity</b>	?90% as determined by SDS-PAGE
<b>Source</b>	This protein contains the mouse TGFR2(Ile24-Asp159) was fused with the C-terminal Fc Tag and expressed in Mammalian cells.
<b>Uniprot ID</b>	<a href="#">Q62312</a>
<b>Expiration Date</b>	6 months from date of receipt.